

SDI

glass ionomer guide

riva



riva self cure

tough, strong, packable glass ionomer restorative

- Self cured, conventional glass ionomer restorative
- Chemically adheres to tooth structure
- Great for minimally invasive dentistry (MID)
- High mechanical strengths – withstands mastication forces
- High fluoride releasing and recharging ability
- Packable
- Suitable for bulk filling
- Fast set initial set time 2.5 min
- Radiopaque
- Non-stick handling

compares to:

Fuji IX GP* Fuji IX GP Extra* Ketac Molar*

Available shades:

A1, A2, A3, A3.5, A4, B2, and B3

Capsules:

- Boxes of 50 capsules in regular and fast set
- Assorted capsule box: 10 each of A1, A2, A3, A3.5 and B2 in regular and fast set

Powder / Liquid:

- Powder (15g) / Liquid (6.9mL/8g) kit
- Powder refill (15g)
- Liquid refill (6.9mL/8g)



riva light cure

the glass ionomer restorative with true esthetics

- Light cured, resin reinforced glass ionomer restorative
- Fabulous esthetics
- Withstands forces of mastication - high compressive strength
- Reduced cracking and failure in flexion – high flexural strengths
- Flexible working time
- Chemically adheres to tooth structure
- Excellent marginal integrity
- Suitable for minimally invasive dentistry (MID)
- Long term wear resistance due to high surface hardness
- Fluoride releasing and recharging
- Varnishing not required

compares to:

Fuji II LC* Photac-Fil Quick*

Available shades:

A1, A2, A3, A3.5, A4, B2, B3, B4, C2, C4 and bleach

Capsules:

- Boxes of 50 capsules
- Assorted capsule box: 10 each of A1, A2, A3, A3.5 and B2

Powder / Liquid:

- Powder (15g) / Liquid (7.2mL/8g) kit
- Powder refill (15g)
- Liquid refill (7.2mL/8g)



riva protect

super high fluoride releasing surface protector and sealant

- ACP nanotechnology – supercharges remineralisation
- Great protector for tooth surfaces, pits and fissures
- Moisture tolerant for use on partially erupted teeth
- Ideal for caries-challenged patients
- Works in a moist environment
- Assists caries prevention – super high fluoride releasing
- Fast and regular set option
- Chemically adheres to tooth structure
- Radiopaque
- Low viscosity for easy placement and pit/fissure penetration

compares to:

Fuji VII* Fuji Triage*

Available shades:

Pink and White

Capsules:

- Boxes of 50 capsules in regular and fast set

Powder / Liquid:

- Powder (15g) / Liquid (9.1mL/10g) kit (includes 10mL Riva Conditioner)
- Powder refill (15g)
- Liquid refill (9.1mL/10g)

Made in Australia by SDI Limited
 Bayswater, Victoria 3153
 Australia 1 800 337 003
 Austria 00800 022 55 734
 Brazil 0800 770 1735
 France 00800 022 55 734
 Germany 0800 100 5759

Ireland 01 886 9577
 Italy 800 780625
 New Zealand 0800 734 034
 Spain 00800 022 55 734
 United Kingdom 00800 022 55 734
 USA & Canada 1 800 228 5166
 www.sdi.com.au

SDI

*Fuji IX GP, Fuji IX GP Extra, Ketac Molar, Fuji II LC, Photac-Fil Quick, Ketac Cem, Fuji I, Fuji Triage, Fuji VII, Miracle Mix, Ketac Silver, Fuji Plus, Relyx Luting Cement are not the registered trademarks of SDI Limited.

SDI

glass ionomer guide



riva silver

mercury free fluoride releasing silver restorative

- Self cured, conventional glass ionomer restorative
- Chemically adheres to teeth
- Condensable - easy handling
- Non-irritating to pulp
- Very high radiopacity
- High wear resistance
- High early strength - perfect for difficult to handle patients
- Great for core build ups
- High fluoride releasing

compares to:

- Miracle Mix*
- Ketac Silver*

Available shade:

Light grey

Capsules:

- Boxes of 50 capsules

Powder / Liquid:

- Powder (32g) / Liquid (8.9mL/10g) kit
- Powder refill (32g)
- Liquid refill (8.9mL/10g)



riva luting

glass ionomer luting cement with fabulous fluoride release

- Self cured, conventional glass ionomer luting cement
- High fluoride releasing
- Chemically adheres to tooth and metal substrates
- Very low film thickness – flows easily between restoration and tooth
- Increased longevity and esthetics - low water solubility

compares to:

- Fuji I*
- Ketac Cem*

Available shade:

Light yellow

Capsules:

- Boxes of 50 capsules

Powder / Liquid:

- Powder (35g) / Liquid (24.3mL/25g) kit
- Powder refill (35g)
- Liquid refill (24.3mL/25g)



riva luting plus

resin modified glass ionomer luting cement

- Self cured, resin modified glass ionomer luting cement
- Fluoride releasing and recharging ability
- Chemically adheres to tooth and substrates
- Withstands forces of mastication - high compressive strength
- Reduced cracking and failure in flexion – high flexural strengths
- Increased longevity and esthetics - low acid water erosion
- Radiopaque

compares to:

- Fuji Plus*
- Relyx Luting Cement*

Available shade:

Light yellow

Capsules:

- Boxes of 50 capsules

Powder / Liquid:

- Powder (25g) / Liquid (8.95mL/10g) kit
- Powder (25g) refills
- Liquid refill (8.95mL/10g)



riva applicator 2

sleek and sophisticated

- Hand contoured design for maximum comfort
- Easy to clean with no sharp edges - autoclavable
- Easy oral access
- Capsule lock provides controlled extrusion
- Superbly engineered stainless steel mechanical parts
- Robust & long lasting

riva coat

- Prevents saliva and water contamination to exposed surfaces of glass ionomer cements during the first 24 hours immediately after placement.
- Bottle refill 5mL (5.75g)

riva conditioner

- Mild solution of polyacrylic acid
- Partially removes the smear layer
- Optimises adhesion of the GIC
- Bottle refill 10mL (10g)

M861002 C



+D036M8610021+

